



## Hallowell EMC Heated Hard Pad for Rodent WorkStand



The Heated Hard Pad for Rodent WorkStand will easily adapt a HEMC Heated Hard Pad to your existing HEMC Rodent WorkStand. Its design will give you the ability to utilize the heated surface while in use on the work stand. The unit installs easily and can be ready for use in minutes with your existing Rodent WorkStand and circulating pump. The unit is also designed so it can be used as a stand-alone Heated Hard Pad. This product is a great versatile addition to your already existing Hallowell EMC products and supplies.

### Installation Instructions

1. Remove existing Alignment Bars on your current WorkStand.
2. Align the Heated Hard Pad up on the top surface of the WorkStand using the slots. The rubber feet on the Hard Pad should come to rest on the WorkStand top surface.
3. Install the 7/8" long Thumbscrews and lock washers (2 each with kit) through the rear side of the WorkStand by inserting them through the slots and threading them into the inserts in the outer border of the Hard Pad. Thread until snug, there's no need to over tighten. This step will secure the Hard Pad to the WorkStand. (See Figure 2)
4. Install the Alignment Bars by using the 1-1/2" long Tee Thumbscrews and flat washers supplied with the kit (2 each with kit). Install Tee Thumbscrews and washers through rear side of WorkStand slots. Thread into Alignment Bars on top surface of the Hard Pad. Leave loosened until ready for use. (See Figure 2)
5. Attach the Hard Pad connectors to the circulating pump. Let the unit heat up and it's ready for use.
6. Reverse these instructions and the Hard Pad and WorkStand can be used as separate units.




Don't forget to check out these related items on our website [www.hallowell.com](http://www.hallowell.com)



**HEMC Heated Hard Pads**

- HEMC Part: 000A2787B  
(21" x 24")
- HEMC Part: 000A2788B  
(15" x 24")
- HEMC Part: 000A2789B  
(9" x 12")
- HEMC Part: 000A4076A  
(6" x 8")



**000A3467 Rodent WorkStand**



**000A3746 Rat Intubation Pack  
000A3747 Mouse Int. Pack**



**000A3748 MDS Otoscope  
w/ Operating Head**



**000A4084 HTP-1500 Pump**