


## Pedestal Stand Attachments Featuring Tool Free Installation Design



Convenient design allows multiple attachment configuration options.  
These items will fit 1-3/4" diameter column pedestals.  
All units can be customized to adapt to a square post.

## **Pedestal Stand Attachments**

**000A6179** – Mayo Instrument Side Tray

**000A6185** – IV Pole Attachment  
(Unit also comes with rolling stand to use as a stand-alone unit)

**000A6186** – E-Cylinder Tank Attachment

**000A6189** – Basket Attachment #1, 10.5” x 5” x 13”

**000A6190** – Basket Attachment #2, 9” x 6” x 11.5”

**000A6191** – Basket Attachment #3, BP Cuff

**000A6192** – O2 DISS Dual Interface Attachment  
(Custom gas interfaces can be made)

**000A6193** – Power Strip Attachment  
(Hospital Grade, 4 outlet w/15’ Cord, 15Amp)

***\*\*Mounting can be customized to other sized columns and square posts\*\****

---

### **No more maintenance calls!**

Our NEW tool free design makes installation a snap with the comfortable knurled hand-tights. The tool free design allows these units to be just as easily removed, adjusted, and relocated by simply undoing the knurled knobs. Units come with easy to read instructions.

