Features Comparison of Junior and Tafonius Large Animal Machines			
Junior sans gas	Junior w/gas	Tafonius (full boat)	
000A5241	000A5242	000A4630	
Ventilator only	Gas machine	Gas machine,	
	and ventilator	ventilator,	
		monitoring and	Features as of 1 JUN 2019
		advanced	
	,	features	
√ 15 seconds	√ 30 seconds	√ 2 seconds	Stainless and Aluminum construction with white powder coated finish Oviet, access to internal planning of the "hallows"
√ V	√ seconds	Z Seconds √	Quick access to internal cleaning of the "bellows" 8" wheels, front wheel locks
	√	√	Accommodates a wide range of patients 100 lbs – 2000 lbs (more or less!)
· √	√	v V	Future, new modes of ventilation as features are software upgradeable
√	√	√	Always powered by two 12v rechargeable batteries, AC keeps them topped off
√	√	√	Servo driven piston, NO bellows
	,	,	• Control Settings for Tidal Volume, Respiratory Rate, Insp-time, MWPL, CPAP/PEEP, Displays the
√ √	√ √	√ √	Resulting MV, I-flow, E-time and I:E ratio
	V	V √	 Displays above parameters when patient is spontaneously breathing; Respiratory monitoring Tidal volumes up to 20 liters in 100ml increments with 25 ml resolution
Future	Future	v V	Delivers compliance compensated TVs insures the patient actually receives the set TV
√ V	√ V	v V	Respiratory Rate 1 – 30 breaths per minute
√	√	√	Maximum Working Pressure Limit adjustable from 10 – 80 cmH20
√	√	√	Implementation of PEEP and CPAP up to 50cmH2O
√ /	√	√	Implementation of Assist Mode Ventilation
√ -/	√ -/	√ -/	Inspiratory flow rates up to 900 lpm Facilitates the implementation of Closed Circuit Anathogia with add a liter feature.
√ √	√ √	√ √	Facilitates the implementation of Closed Circuit Anesthesia with add-a-liter feature Scavenging flowmeter
		V √	Backup scavenging system
	√	√	Uses NO driving gas to compress the "bellows"; a major ongoing cost savings
√	√	√	Oxygen is used only for patient uptake
Optional	√	√	Airway pressure gauge
Additional	Additional	Additional	• IV pole
	√ -/	√ -/	22mm Foal Breathing Circuit included Lightweight Light free gloop breathing tubes
	√ √	√ √	Lightweight, kink free clear breathing tubes Accommodates both funnel fit and Bivona/Drager style ET tube adapters
			Absorber pan holds 6.5 liters, 11.7 lbs, 5.3 kg of absorbent
	√	v V	Oxygen flowmeter for manual use
	√	√	Oxygen flush button for manual use
	Additional	√	Hospital line pressure gauge
	√	√	Low oxygen supply pressure alarm
	Additional	Additional	Optional N2O and/or Air flowmeters
	n*	√ √	Bronchodilator injection port Pullout swiveling keyboard tray
	n		Dual mount Select-a-Tec™ vaporizer back bar
	n	v V	Swivel valve block permits the breathing tubes to exit to the right or the left
	n	√	Stainless hoop handle and bumpers
	n	√	GFI protected with 3 auxiliary AC power outlets
	n	√ ,	Panel mount e-cylinder pressure gauge
	n n	√ √	"Crash cart" drug drawer below work surface Absorbent can be changed in seconds
	n	V √	Two work lights and a built in maintenance tool kit
	n		Accommodates two oxygen e-cylinders for backup, temporary disconnect or field use
Additional	Additional	Included	Optional oxygen e-cylinder mount for backup, temporary disconnect or field use
		√	• Footprint 30"wide x 36"deep, height is 72"
		√	Redundant computer systems
Main Cantual	Main Control	Main Control	15" 1028x768 touchscreen computer w/GUI interface and 160 Gigabyte HD
Main Control	Main Control	V Sackup Control	Microprocessor with TV, RR, Itime, MWPL, CPAP/PEEP and Assist mode controls Settable In inspiratory pages 0.40% of the Itime
			Settable Ip, inspiratory pause, 0-40% of the Itime Settable FiO2, 30% – 100 % (requires Air Option above)
		v V	Choice of automatic or manual modes for FiO2 and anesthetic delivery
		√	• Fully integrated monitoring of ECG, IBP, SpO2, Temp, Insp/Expired CO2, O2, Agent, N2O
		√	Sidestream Gas Analyzer standard
		√	Paramagnetic oxygen sensor w/ above, no galvanic cells to worry about or replace!
		√	High and low alarms for all parameters
		√ √	Anesthesia record printout for clinical use Data collection to Excell ^M for recognish
-		V Future	Data collection to Excel™ for research User programmability for the researcher
		ruture √	Unlimited presets for recalling all settings and alarms for different doctors or types of cases
			Trend display of all measured parameters
		v V	User configurable display of traces and trends
		√	Freeze traces function
		√	Alarm silence/locked/disabled function
√	√		Alarm silence function only
	- /	l	• Footprint 24"wide x 28"deep, height is 50"
√	√ ^ dditional		
√ Additional	Additional		ZOO PACKAGE includes collapsable lifting arms & mayo trays for ease of use in the field
√ Additional Optional			ZOO PACKAGE includes collapsable lifting arms & mayo trays for ease of use in the field VML by MidMark (Matrx) gas machine head and mounting pole
√ Additional		Depends on	ZOO PACKAGE includes collapsable lifting arms & mayo trays for ease of use in the field

^{*} Note: It is possible to upgrade Junior with gas to have many of the feature of Tafonius. The features that would not be available are indicated by the symbol "n" in the table above. Pricing for features indicated as "Optional" has not been determined. DOCA4958A