



HALLOWELL EMC

000A6558

MRI Ventilator Conversion Kit w/O2 DISS Inline Pressure Gauge



Introducing the new Hallowell MRI Conversion Kit.

This unit easily adapts to your existing HEMC Pro Series Model ventilators. Along with an updated design, the kit comes with 25' of extension tubing to reach your MRI environment, and the side mount panel has an O2 pressure gauge that will allow the input supply gas to be monitored at the unit.



This kit is designed for easy installation to your existing ventilator. Simply remove the left side instruction panel from the unit. Using the 4 screws removed, mount the panel gauge.



Full hook-up instructions available with purchased unit or in our Document Library at www.hallowell.com



Comes with optional use rubber feet for base bottom. Useful if the unit is not secured to an anesthesia machine.