



Specifications

Battery	One AA battery, lasts approx 360 hours
Weight	108 g (3.8 oz)
Size	127 x 62 x 28mm (5 x 2.4 x 1.1in)
Measurement Units	Flow rate (mL/h), drops per minute (dp/m) and total volume (mL)
Alert	75 dB at 10cm approx. Sounds at $\pm 13\%$ rate change or when flow stops
Minimum Chamber Diameter	13 mm (0.51in)
Maximum Chamber Diameter	21.1 mm (0.83in)
Compatible Drip Sets	10, 15, 20, 60 gtt/mL
Accuracy	$\pm 1\%$ drop detection
Warranty	1 year limited warranty

Standards Compliance

Radiated Emissions	CISPR 11:2015 +A1:2016
Ingress Protection	IP22
Electromagnetic compatibility	IEC 61000-4
Regulatory	ISO 13485:2016

Additional Warnings for Customers in the USA

The following warnings apply only for customers intending to use the DripAssist device within the United States of America:

- Only trained, qualified personnel should operate the DripAssist.
- The DripAssist is not intended for home use.
- Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.