

# HALLOWELL EMC

## 000A6770 - Veterinary Mobile Anesthesia Station w/Storage

Safely maneuver around a facility with a complete and secure workstation! The stable cart design takes up less floor space than a rolling pedestal stand and allows a vaporizer, ventilator and/or a monitoring system to safely travel. Locking casters are included.

The mobile cart is the perfect height, with both a pullout drawer and lower storage cabinet that lock. There is the convenience of a fold over hinged shelf for placement of a ventilator or monitoring system. The workstation surface also has extra space to work on.

The unit is equipped with a single station Select-A-Tec vaporizer mount that can easily be changed out with the cagemount adaptor option available. The mobile cart comes fully assembled to your facility. Add your optional components and you're good to go!



Customizable to **YOUR** requirements

Ships Fully Assembled

For pricing and more information  
please contact us at  
[sales@hallowell.com](mailto:sales@hallowell.com)

## **000A6770 - Veterinary Mobile Anesthesia Station w/Storage**

### **Standard Features:**

- 0-10 LPM Parker O<sub>2</sub> Flowmeter
- Securely placed Flush Valve (calibrated at 40 LPM)
- 1500 ml CO<sub>2</sub> absorber designed for easy and fast replenishing of sodalime. Condensation trap to prevent water build-up in sodalime. Easy and efficient gas sealing. The absorber hosts a sampling port for both the Expiration and Inspiration ports. The manometer is directly ported to the inhale side.
- Leak test port on the absorber, simply plug the wye connection into the port while leak testing the system.
- Convenient Fold Over Hinged Shelf
- Secure and easy to change tubing connections.
- Select-A-Tec single station manifold.
- Sturdy assembled cart w/locking casters.
- O<sub>2</sub> DISS supply line inlet.
- Pressure gauge to monitor supply line input.
- Additional oxygen port for easy hookup of ventilator.
- Sleek power bar design that includes 4 USB A ports.

### **Accessories Supplied:**

- 4' Sampling lines (2)
- Breathing Circuit
- f/Air canister holder w/square post clamp
- f/Air canister
- 19mm Purple Scavenger Tube (5')
- 4' O<sub>2</sub> DISS Supply Hose, Grn.
- O<sub>2</sub> DISS Male Coupler (Everyone should have one of these handy)
- 2 Liter EVAC Purple Reservoir Bag
- Set of 4 Thumbscrews to secure ventilator to fold-over shelf.

### **Optional Items Available:**

- 2024Pro Ventilator
- TEC7 ISO Select-A-Tec mount, Quick Fill Vaporizer
- Cagemount Vaporizer Adapter Kit
- Vac/Scav Interface
- Basket attachments
- IV Pole attachment

