

Model 2002IE Pro Series

Front & Rear



Model 2024IE Pro Series

Front & Rear



Summary of upgrades between the 2002 Pro and 2024 Pro Series.

The primary function of the newest unit provides the end user more precision digital ventilation. The 2024 Pro Series is a microprocessor-controlled unit providing the latest in electronic updates. The function of this version remains primarily the same as the previous Pro Series.

The following is a quick reference list of additions/upgrades between the 2002 Pro Series and 2024 Pro Series, all of these can be found in the **User Manual** for complete explanation and proper usage.

Overall unit:

The chassis is now powder coated white with black silkscreened lettering. Easier to clean with no overlay to damage. The covers are still powdered blue. Overall better medical look.

Front of unit:

- The Power Switch is now a backlit LED pushbutton switch.
- Added function of a Standby/Ventilate Switch w/LED.
- Rate Control w/digital display.
- Volume Control valve upgrade.
- MWPL Control w/precision digital display.
- Inspiratory Hold pushbutton with LED aid.
- Updated Alarm LEDs added to aid with unit operation. (MWPL, Low Breathing & Low Gas)

Rear of unit:

- Medical safety approved 12 volts, 2.1mm plug-in power supply brick.
- 2amp Fuse Holder
- Optional feature ports for Manual Remote Inspiratory Hold and Oxygen Flow Interrupter.

All previous consumables fit the 2024 Pro Series units.